

# Piezo Mini

## MINI DENTAL PIEZO SURGERY

Minimize and lightweight  
dual-screen design for smart dual control



# MINI DESIGN

## COMPACT & EFFICIENT

Balancing minimally invasive techniques & power output

# 5"

MAIN SCREEN



Full-touch interface for intuitive operation

# 1.6"

ROUND SCREEN



Real-time data display and adjust

### OPTIMIZED DIMENSIONS

perfectly designed for a streamlined workspace

### SLIM HANDPIECE

Operates at full power for 10 mins with a handpiece temperature of only 33°C, ensuring safety and comfort



Bone Mode



Endo Mode



Perio Mode



Custom Mode



Clean Mode

### PARAMETER

**VOLTAGE** 100~240V 50/60Hz 120VA

**RIDAL AMPLITUDE** 20-200µm

**FUSES** 2×F10A 250V

**LATERAL AMPLITUDE** ≤30µm

**FREQUENCY** 26~38Hz

**WATER FLOW** 60-100ml

<https://www.alandental.com>

PLUWA PLUS POWER SG60P  
SN: S0524JUA11823